

PEBAX® 7233 SA 01 MED

Polyether block amide Pebax® 7233 SA 01 MED is a thermoplastic elastomer made of flexible polyether and rigid polyamide. This grade offers the highest quality and it is specially designed to meet the stringent requirements of the medical applications such as minimally invasive devices. Pebax® 7233 SA 01 MED also offers an excellent combination of properties such as: kink resistance, low friction coefficient and superior dynamic response. Upon request, letters regarding USP Class VI compliance can be provided.

MAIN CHARACTERISTICS

Property	Typical Value	Unit	Test Method
Density	1.01	g/cm ³	ISO 1183
Water Absorption at Equilibrium At 20°C and 50 % R.H.	0.7	%	ISO 62
Water Absorption At 23°C and 24 h in water	0.9	%	
Melting Point	174	°C	ISO 11357
Vicat Point Under 1 daN	164	°C	ISO 306
Shrinkage (after 24 h, 4 mm, mold at 40°C) // ⊥	1.2 1.5	% %	Internal method
Hardness (*) Instantaneous After 15 s	69 61	Shore D Shore D	
Tensile Test (*) Stress at Break Strain at Break	56 >300	MPa %	ISO 527
Flexural Modulus (*)	510	MPa	ISO 178
Charpy Impact (*) Unnotched 23°C Unnotched -30°C V-notched 23°C V-notched -30°C	No break No break 15 (c) 10 (c)	kJ/m² kJ/m² kJ/m² kJ/m²	ISO 179

(*) Samples conditioned 15 days at 23°C - 50 % R.H.

(c) = Complete Break

PEBAX® 7233 SA 01 MED

MAIN APPLICATIONS

- Tubings like angiography and angioplasty catheters.
- Flexible injected parts.

PROCESSING CONDITIONS

Conditions	Typical values
Extrusion Melt Temperature (Min / Recommended / Max)	220°C / 235°C / 250°C
Injection Melt Temperature (Min / Recommended / Max)	230°C / 260°C / 290°C
Mold Temperature	25 – 60°C
Drying (only necessary for bags opened for more than two hours) Time Temperature	5 - 7 hours 70 - 80°C

PACKAGING

This grade is delivered dried in sealed packaging (20 kg bags) ready to be processed.