+135-3858-6433 (GuangDong) +188-1699-6168 (ShangHai) +852-6957-5415 (HongKong)

We Connect Science



LUPOY PC2001-20

Injection Molding & Extrusion Grade, PC

Description

Application

FDA compliance¹⁾, Good mold release, Low viscosity, Bio-compatibility ISO 10993 certified EtO and steam sterilization at 121°C

Medical Devices (Syringes, Dialyzers, Connectors)

Ver. 23.07

Properties	Test Condition	Test Method	Unit	Typical Value
Physical				
Specific Gravity		ISO 1183	g/cm ³	1.20
Molding Shrinkage		ISO 294-4	mm/mm	0.005~0.007
Melt Flow Rate	300℃/1.2kg	ISO 1133	cm ³ /10min	19
Optical				
Refractive Index, np		ISO 489		1.586
Light Transmittance		ASTM D 1003	%	89
Haze		ASTM D 1003	%	<1.5
Tensile Strength (@ Yield) Tensile Flongation (@ Break)	50mm/min 50mm/min	ISO 527 ISO 527	MPa %	60 >100
Mechanical				
Tensile Elongation (@ Break)	50mm/min	ISO 527	%	>100
Flexural Strength	2mm/min	ISO 178	MPa	96
Flexural Modulus	2mm/min	ISO 178	MPa	2,300
Charpy Impact Strength		ISO 179		
(Notched)	23℃		KJ/m ²	65
Izod Impact Strength		ISO 180		
(Notched)	23℃		KJ/m ²	65
Thermal				
Heat Deflection Temperature		ISO 75		
(unannealed, Flatwise)	0.45MPa		${\mathbb C}$	131
(unannealed, Flatwise)	1.8MPa		$^{\circ}$	121
Vicat Softening Point	5kg, 50°C/h	ISO 306	°C	142

comply with the U.S. Food, Drug, and Cosmetic Act and Food Additive Regulations 21 CFR 177.1580 and E.U. Food Contact Regulations

1) When used unmodified for the manufacture of food contact articles LUPOY PC2001-20 Polycarbonate resins

Note) Typical values are only for material selection purpose, and variation within normal tolerances are for various colors.

Values given should not be interpreted as specification and not be used for part or tool design.

All properties, except melt flow rate are measured on injection molulded specimens and after 48 hours storage at 23 °C, 50% relative humidty.

We Connect Science



LUPOY PC2001-20

Injection Molding & Extrusion Grade, PC

Description

Application

FDA compliance¹⁾, Good mold release, Medical Devices

Low viscosity, Bio-compatibility ISO 10993 certified (Syringes, Dialyzers, Connectors)

EtO and steam sterilization at 121°C

Ver. 23.07

Electrical

GWT 2.0 mm, 5 second		IEC 60695-2-13	${\mathbb C}$	850
Comparative Tracking Index(CTI)	Solution A	IEC 60112	Volts	250
Volume Resistivity	23°C	ASTM D257	Ohm-m	2x10 ¹⁷
Dielectric Strength	23℃	ASTM D149	KV/mm	17
Dielectric Constant (60 Hz)	23°C	ASTM D150		3
Dissipation Factor (60 Hz)	23℃	ASTM D150		0.001
Comparative Tracking Index (2.00 mm)	23℃	IEC 60112	V	250

Note) Typical values are only for material selection purpose, and variation within normal tolerances are for various colors.

Values given should not be interpreted as specification and not be used for part or tool design.

All properties, except melt flow rate are measured on injection molulded specimens and after 48 hours storage at 23 °C, 50% relative humidty.

Processing Guide (Injection Molding)

Process	ing Parameters	Unit	Value
Drying Temperature		°C	120
Drying Time		hrs	4
Maximum Moisture Content		%	0.02
Melt Temperature		$^{\circ}\! \mathbb{C}$	280 ~ 300
Cylinder Temperature	Rear	$^{\circ}$	265 ~ 285
	Middle	$^{\circ}$ C	270 ~ 290
		$^{\circ}$	275 ~ 295
Nozzle Temperature		$^{\circ}$	280 ~ 300
Mold Temperature		°C	80 ~ 120
Back Pressure	_	kg/cm ²	10 ~ 40
Screw Speed	_	rpm	40 ~ 70

Note) Back Pressure & Screw Speed are only mentioned as general guidelines.

These may not apply or need adjustment in specific situations such as low shot sizes, thin wall molding and gas-assist molding.